




REV.	ECO NO.	DESCRIPTION
A	000023	RELEASE TO PRODUCTION.



NOTE:

MAX VOLUME OF RESERVOIR 90 ml

**FOR REFERENCE ONLY**

**INNOVATIVE**  
CUSTOM MICROPLATE  
MICROPLATES

1998 WESTOVER RD.  
CHICOPEE, MA 01022  
TEL. (413) 593-5300  
FAX (413) 593-5900

© 2001  
CONFIDENTIAL  
TRADE SECRETS  
PATENTS ISSUED  
AND PENDING  
DO NOT COPY  
OR DISCLOSE  
WITHOUT WRITTEN  
PERMISSION OF  
INNOVATIVE MICROPLATE

TOLERANCES UNLESS  
OTHERWISE SPECIFIED  
XX ± .03  
XXX ± .015  
XXXX ± .0005  
FRACTION: ± 1/64  
ANGULAR: ± 0°15'  
RADIUS 0.015" UNLESS  
OTHERWISE SPECIFIED  
BREAK ALL SHARP EDGES

MATERIAL: RESPPCP  
NATURAL POLYPROPYLENE  
CONCENTRATE:

DATE: 02/12/02  
DRAWN BY: YELENA SUKONIK

PRT NAME:

RESERVOIR, LOW PROFILE, FLAT BOTTOM, NATURAL PP

CATALOG NO.: S30023

APPR BY:

PRT SIZE:  
A

SCALE:  
N15

SHEET

1 OF 1

REV.  
A