




REV.	ECO NO.	DESCRIPTION
A	000023	RELEASE TO PRODUCTION.



- NOTES:
1. MAX VOLUME OF RESERVOIR 282 ml
 2. MAX VOLUME OF WELL 0.007 ml



FOR REFERENCE ONLY

INNOVATIVE
CUSTOM MICROPLATE
MICROPLATES

1998 WESTOVER RD.
CHICOPEE, MA 01022
TEL. (413) 593-5300
FAX (413) 593-5900

© 2001 CONFIDENTIAL TRADE SECRETS PATENTS ISSUED AND PENDING DO NOT COPY OR DISCLOSE WITHOUT WRITTEN PERMISSION OF INNOVATIVE MICROPLATE	TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .0005 FRACTION: ± 1/64 ANGULAR: ± 0°15' RADIUS 0.015" UNLESS OTHERWISE SPECIFIED	MATERIAL: RESPPCP NATURAL POLYPROPYLENE CONCENTRATE:
	DATE: 02/12/02 DRAWN BY: YELENA SUKONIK	PRT NAME: 384 WELL, RESERVOIR, PYRAMID BOTTOM, NATURAL PP

APPR BY:	PRT SIZE:	SCALE:	SHEET	OF	REV.
	A	NIS	1	1	A