


REV.	ECO NO.	DESCRIPTION
A	000377	RELEASED FOR PRODUCTION



FOR REFERENCE ONLY

INNOVATIVE
CUSTOM **MICROPLATE**
MICROPLATES

1998 WESTOVER RD.
CHICOPEE, MA 01022
TEL. (413) 593-5300
FAX (413) 593-5900

© 2001
CONFIDENTIAL
TRADE SECRETS
PATENTS ISSUED
AND PENDING
DO NOT COPY
OR DISCLOSE
WITHOUT WRITTEN
PERMISSION OF
INNOVATIVE MICROPLATE

TOLERANCES UNLESS
OTHERWISE SPECIFIED
XX ± .03
XXX ± .015
XXXX ± .0005
FRACTION: ± 1/64
ANGULAR: ±0°15'
RADIUS 0.015" UNLESS
OTHERWISE SPECIFIED
BREAK ALL SHARP EDGES

MATERIAL: RESPPCP
NATURAL POLYPROPYLENE
CONCENTRATE:

DATE: 7/18/05
DRAWN BY: KN

PRT NAME:
384 WELL, PARTITIONED RESERVOIR, LOW PROFILE
PYRAMID BOTTOM, W/WAVE BEAKERS, NATURAL PP

CATALOG NO.: S30047

APPR BY:	PRT SIZE: A	SCALE: NTS	SHEET 1 of 1	REV. A
-----------------	-----------------------	----------------------	---------------------	------------------