


REV.	ECO NO.	DESCRIPTION
A	000023	RELEASE TO PRODUCTION.



- NOTES:
1. MAX VOLUME OF RESERVOIR 92 ml
 2. MAX VOLUME OF WELL 0.007 ml

FOR REFERENCE ONLY

INNOVATIVE
CUSTOM MICROPLATE
MICROPLATES

1998 WESTOVER RD.
CHICOPEE, MA 01022
TEL. (413) 593-5300
FAX (413) 593-5900

MATERIAL: RESPPCP
NATURAL POLYPROPYLENE
CONCENTRATE:

TOLERANCES UNLESS OTHERWISE SPECIFIED
XX ± .03
XXX ± .015
XXXX ± .0005
FRACTION: ± 1/64
ANGULAR: ± 0°15'
RADIUS 0.015" UNLESS OTHERWISE SPECIFIED
INNOVATIVE MICROPLATE BREAK ALL SHARP EDGES

DATE: 02/12/02
DRAWN BY: YELENA SUKONIK

PRT NAME:
384 WELL, RESERVOIR, LOW PROFILE, PYRAMID BOTTOM, NATURAL PP
CATALOG NO.: S30022

APPR BY:
PRT SIZE: A
SCALE: NTS

SHEET 1 **OF** 1
REV. A