




REV.	ECO NO.	DESCRIPTION
A	000022	RELEASE TO PRODUCTION.



NOTES:  
 1. MAX VOLUME OF RESERVOIR 86 ml  
 2. MAX VOLUME OF WELL 0.064 ml



<b>FOR REFERENCE ONLY</b> <b>INNOVATIVE</b> CUSTOM MICROPLATE MICROPLATES 1998 WESTOVER RD. CHICOPEE, MA 01022 TEL. (413) 593-5300 FAX (413) 593-5900		TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .0005 FRACTION: ± 1/64 ANGULAR: ± 0°15' RADIUS 0.015" UNLESS OTHERWISE SPECIFIED INNOVATIVE MICROPLATE BREAK ALL SHARP EDGES	<b>MATERIAL:</b> RESPPCP NATURAL POLYPROPYLENE <b>CONCENTRATE:</b>	<b>PRT NAME:</b> 96 WELL, RESERVOIR, LOW PROFILE, PYRAMID BOTTOM, NATURAL PP	<b>SCALE:</b> NTS	<b>SHEET</b> 1 OF 1	<b>REV.</b> A
		<b>DATE:</b> 01/22/02 <b>DRAWN BY:</b> YELENA SUKONIK	<b>CATALOG NO.:</b> S30018	<b>PRT SIZE:</b> A	<b>APPR BY:</b>	<b>APPR BY:</b>	<b>SHEET</b> 1 OF 1