




REV. A	ECO NO. 000020	DESCRIPTION RELEASE TO PRODUCTION.
--------	----------------	---------------------------------------



SEE NOTE#2



- NOTES:
 1. MAX VOLUME OF RESERVOIR 290 ml
 2. MAX VOLUME OF ROW 0.825 ml

INNOVATIVE CUSTOM MICROPLATE MICROPLATES 1998 WESTOVER RD. CHICOPEE, MA 01022 TEL. (413) 593-5300 FAX (413) 593-5900	© 2001 CONFIDENTIAL TRADE SECRETS PATENTS ISSUED AND PENDING DO NOT COPY OR DISCLOSE WITHOUT WRITTEN PERMISSION OF INNOVATIVE MICROPLATE	TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .0005 FRACTION: ± 1/64 ANGULAR: ± 0°15' RADIUS 0.015" UNLESS OTHERWISE SPECIFIED	MATERIAL: RESPPCP NATURAL POLYPROPYLENE CONCENTRATE:		PRT NAME: 8 ROW, RESERVOIR, PYRAMID BOTTOM, NATURAL PP	REV. A
			DATE: 12/14/01 DRAWN BY: YELENA SUKONIK	PRT NO.: S30017	SCALE: NTS	SHEET 1 OF 1